

## Appendix D



Directions with angle diagram assists guests to learn how to align and read the goniometer. The goniometer has the stationary arm color coded and labeled. Movable arm and axis is labeled as well. This goniometer is taped so the arms do not move and can be placed on top of the diagram. See following page for example of how the directions were written.



Additional diagrams with arrows demonstrating the direction of the angle are placed next to the above activity. This allows the guest to practice the skill. The goniometer is labeled and the arms are movable. The degrees of the angle are written on the back of the diagrams so the guests can compare their readings with the correct answer.